

Product Datasheet

ACTH antibody (orb302731)

Description	Mouse monoclonal antibody to ACTH
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide
Concentration	1.0mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.