

Product Datasheet

mast cell tryptase antibody (orb302748)

Description	Mouse monoclonal antibody to mast cell tryptase
Species/Host	Mouse
Reactivity	Canine, Feline, Human, Monkey, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide
Concentration	1.0 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,,,
Expiration Date	12 months from date of receipt.