

Product Datasheet

TGF beta antibody (orb302795)

Description	Mouse monoclonal antibody to TGF beta
Species/Host	Mouse
Reactivity	Human, Mammal, Mink, Mouse, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	PBS
Concentration	1.0 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	
Expiration Date	12 months from date of receipt.