

Product Datasheet

CD326 antibody (orb302797)

Description	Mouse monoclonal antibody to CD326
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide and 0.5% rAlbumin
Concentration	0.1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.