

## Product Datasheet

### macrophages antibody (orb302798)

|                            |                                                                                                                            |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to macrophages                                                                                   |
| <b>Species/Host</b>        | Mouse                                                                                                                      |
| <b>Reactivity</b>          | Bovine, Canine, Deer, Equine, Feline, Goat, Guinea pig, Hippopotamus, Human, Marmoset, Mink, Porcine, Primate, Rabbit, Rat |
| <b>Conjugation</b>         | Unconjugated                                                                                                               |
| <b>Tested Applications</b> | FC, IHC-Fr, IHC-P                                                                                                          |
| <b>Immunogen</b>           | Human monocytes.                                                                                                           |
| <b>Preservatives</b>       | PBS with 0.09% sodium azide                                                                                                |
| <b>Concentration</b>       | 1.0 mg/ml                                                                                                                  |
| <b>Storage</b>             | 2-8°C                                                                                                                      |
| <b>Note</b>                | For research use only                                                                                                      |
| <b>Isotype</b>             | IgG1                                                                                                                       |
| <b>Clonality</b>           | Monoclonal                                                                                                                 |
| <b>Dilution Range</b>      | „                                                                                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                            |