

Product Datasheet

CD63 antibody (FITC) (orb302804)

Description	Mouse monoclonal antibody against CD63 conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	Please contact us for more information.
Preservatives	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.