

Product Datasheet

LRP antibody (orb302828)

Description	Mouse monoclonal antibody to LRP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Preservatives	supernatant with 1% rAlbumin and 0.1% sodium azide
Concentration	250 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.