

Product Datasheet

P193/VPARP antibody (orb302841)

Description	Mouse monoclonal antibody to P193/VPARP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 1% rAlbumin and 0.1% sodium azide
Concentration	250 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG2b k
Clonality	Monoclonal
Dilution Range	IHC-P (1:10-1:20),,,
Expiration Date	12 months from date of receipt.