

Product Datasheet

MRP14 antibody (orb302843)

Description	Mouse monoclonal antibody to MRP14
Species/Host	Mouse
Reactivity	Bovine, Canine, Equine, Feline, Goat, Guinea pig, Human, Monkey, Porcine, Primate, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 ug/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	
Expiration Date	12 months from date of receipt.