

Product Datasheet

MVP antibody (orb302847)

Description	Mouse monoclonal antibody to MVP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P, IP, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 ug/ml
Storage	2-8°C
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Dilution Range	,,,,
Expiration Date	12 months from date of receipt.