

Product Datasheet

MRP6 antibody (orb302849)

Description	Mouse monoclonal antibody to MRP6
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, WB
Immunogen	Please contact us for more information.
Preservatives	supernatant with 1% rAlbumin and 0.1% sodium azide
Concentration	100 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,IHC-fr (1:20-1:50),,WB (1:20-1:50)
Expiration Date	12 months from date of receipt.