

Product Datasheet

Breast Cancer Resistance Protein antibody (orb302851)

Description	Mouse monoclonal antibody to Breast Cancer Resistance Protein
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, WB
Immunogen	Please contact us for more information.
Preservatives	cell culture supernatant with 0.7% rAlbumin and 0.1% sodium azide
Concentration	250 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	ICC (1:20-1:50),IHC-fr (1:20),WB (1:50),
Expiration Date	12 months from date of receipt.