

Product Datasheet

CD120a (phospho-55/60) antibody (orb302856)

Description	Mouse monoclonal antibody to CD120a (phospho-55/60)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FC, IHC-Fr, IP, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.1% rAlbumin and 0.02% sodium azide
Concentration	100 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	,,,,,
Expiration Date	12 months from date of receipt.