

Product Datasheet

CD120b (phospho-75/80) antibody (orb302858)

Description	Mouse monoclonal antibody to CD120b (phospho-75/80)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FC, IF, IP
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.1% rAlbumin and 0.02% sodium azide
Concentration	100 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,,,,
Expiration Date	12 months from date of receipt.