

Product Datasheet

hRrp4 Immunoglobulin antibody (orb302878)

Description	Mouse monoclonal antibody to hRrp4 Immunoglobulin
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Please contact us for more information.
Preservatives	lyophilized in a 10 mM ammonium bicarbonate buffer. Each vial contains 2 mg rAlbumin
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.