

Product Datasheet

Citrullinated Fibrinogen antibody (orb302887)

Description	Mouse monoclonal antibody to Citrullinated Fibrinogen
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Please contact us for more information.
Preservatives	lyophilized in a 10 mM ammonium bicarbonate buffer. Each vial contains 2 mg rAlbumin
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	IHC(1:100-1:500),WB (1:200-1:1000),ELISA(1:2000-1:4000)
Expiration Date	12 months from date of receipt.