

Product Datasheet

Mucin-1 antibody (orb302898)

Description	Mouse monoclonal antibody to Mucin-1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	Contains 0.09% sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,IHC-P (1:20),
Expiration Date	12 months from date of receipt.