

Product Datasheet

S100A8 antibody (orb302904)

Description	Mouse monoclonal antibody to S100A8
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC-Fr, WB
Immunogen	Please contact us for more information.
Preservatives	cell culture supernatant with 0.7% rAlbumin and 0.1% sodium azide
Concentration	250 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	,IHC-fr (1:20),WB (1:50)
Expiration Date	12 months from date of receipt.