

Product Datasheet

Langerh. cells antibody (orb302911)

Description	Mouse monoclonal antibody to Langerh. cells
Species/Host	Mouse
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	cell culture supernatant with 0.7% rAlbumin and 0.1% sodium azide
Concentration	500 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
Dilution Range	IHC-fr (1:100-1:1000),
Expiration Date	12 months from date of receipt.