

Product Datasheet

Interleukin-10 antibody (orb302925)

Description	Rabbit polyclonal antibody to Interleukin-10
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA, IF, IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.1% rAlbumin and 0.02% sodium azide
Concentration	100 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	Ig
Clonality	Polyclonal
Dilution Range	,,,,,
Expiration Date	12 months from date of receipt.