

Product Datasheet

Interferon gamma antibody (orb302928)

Description	Rabbit polyclonal antibody to Interferon gamma
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA, IHC-Fr, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.1% rAlbumin and 0.02% sodium azide
Concentration	100 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	Ig
Clonality	Polyclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.