

Product Datasheet

Laminin antibody (orb302941)

Description	Rabbit polyclonal antibody to Laminin
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	Contains 0.09% sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Clonality	Polyclonal
Dilution Range	,IHC-P (1:50 - 1:100),
Expiration Date	12 months from date of receipt.