

Product Datasheet

HIST1H2BG antibody (orb30299)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	H2BC4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	13906
Uniprot ID	P62807
NCBI	NP_003514.2 , NP_003513.1 , NP_066407.1 , NP_003509.1 , NP_003517.2 , NP_003516.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

