

Product Datasheet

VMAT2 antibody (orb303049)

Description	Rabbit polyclonal antibody to VMAT2
Species/Host	Rabbit
Reactivity	Canine, Guinea pig, Hamster, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	Lyophilized; reconstitute in 100 µl dist. water
Concentration	n/a
Storage	2-8°C
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-fr (1:800-1:1,
Expiration Date	12 months from date of receipt.