

Product Datasheet

Neuropeptide Y 1-20 antibody (orb303064)

Description	Rabbit Polyclonal antibody to Neuropeptide Y 1-20
Species/Host	Rabbit
Reactivity	Hamster, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	Lyophilized; reconstitute in 100 µl dist. water
Concentration	n/a
Storage	2-8°C
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-fr (1:200-1:400)
Expiration Date	12 months from date of receipt.