

Product Datasheet

VGAT antibody (orb303070)

Description	Rabbit polyclonal antibody to VGAT
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	Lyophilized; reconstitute in 100 µl dist. water
Concentration	n/a
Storage	2-8°C
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-fr (1:1,,
Expiration Date	12 months from date of receipt.