

Product Datasheet

Helospectin antibody (orb303071)

Description	Rabbit polyclonal antibody to Helospectin
Species/Host	Rabbit
Reactivity	Canine, Feline, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	Lyophilized; reconstitute in 100 µl dist. water
Concentration	n/a
Storage	2-8°C
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-fr (1:400-1:800)
Expiration Date	12 months from date of receipt.