

Product Datasheet

Pancreastatin antibody (orb303073)

Description	Rabbit polyclonal antibody to Pancreastatin
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Preservatives	Lyophilized; reconstitute in 100 µl dist. water
Concentration	n/a
Storage	2-8°C
Note	For research use only
Clonality	Polyclonal
Dilution Range	,IHC-P (1:200-1:1,
Expiration Date	12 months from date of receipt.