

Product Datasheet

CGRP antibody (orb303085)

Description	Rabbit polyclonal antibody to CGRP
Species/Host	Guinea pig
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Preservatives	Lyophilized; reconstitute in 100 µl dist. water
Concentration	n/a
Storage	2-8°C
Note	For research use only
Clonality	Polyclonal
Dilution Range	,IHC-P (1:800-1:1,,
Expiration Date	12 months from date of receipt.