

Product Datasheet

Adrenomedullin antibody (orb303086)

Description	Rabbit polyclonal antibody to Adrenomedullin
Species/Host	Guinea pig
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	Lyophilized; reconstitute in 100 µl dist. water
Concentration	n/a
Storage	2-8°C
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-fr (1:400 - 1:800)
Expiration Date	12 months from date of receipt.