

## Product Datasheet

### HCG antibody (orb303132)

|                            |                                                                       |
|----------------------------|-----------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to HCG                                     |
| <b>Species/Host</b>        | Rabbit                                                                |
| <b>Reactivity</b>          | Human                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                          |
| <b>Tested Applications</b> | IHC-P                                                                 |
| <b>Immunogen</b>           | Please contact us for more information.                               |
| <b>Preservatives</b>       | Tris Buffer, pH 7.3-7.7, containing 1% rAlbumin and 0.1% Sodium Azide |
| <b>Concentration</b>       | n/a                                                                   |
| <b>Storage</b>             | 2-8°C                                                                 |
| <b>Note</b>                | For research use only                                                 |
| <b>Clonality</b>           | Polyclonal                                                            |
| <b>Dilution Range</b>      | IHC-P (1:50-1:250)                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                       |