

Product Datasheet

**FVIIa | Coagulation factor VIIa (500 µg)
(orb303161)**

Description

FVIIa | Coagulation factor VIIa (500 µg)

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant human Factor VIIa (NovoSeven , Novo Nordisk A/S Denmark)
Form/Appearance	Lyophilized
Storage	Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them, to avoid any losses that might occur from material adhering to the cap or sides of the tube.
Note	For research use only
Application notes	This antibody is purified using Protein G, to obtain a total IgG fraction
Clonality	Polyclonal
Purity	Immunogen affinity purified IgG in PBS pH 7.4.
MW	52 kDa (this is MW of Factor VIIa)
Uniprot ID	Q9SXJ7 , Q9FI56
Dilution Range	I-ELISA: 1: 15 000, IL: 14 µg/ml
Expiration Date	12 months from date of receipt.