

Product Datasheet

NCAM 2 antibody (orb303188)

Description

Rabbit polyclonal antibody to NCAM 2

Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant polypeptide corresponding to amino acid residues 442-685 of mouse N-CAM O35136
Form/Appearance	Lyophilized
Storage	Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.
Note	For research use only
Application notes	N-CAM 2 is a cell adhesion molecule belonging to the immunoglobulin superfamily. It is differentially expressed in the developing and mature olfactory epithelium and is involved in the guidance of developing olfactory sensory neurons to the glomeruli of the olfactory bulb. Alternative names: RNCAM, OCAM.
Clonality	Polyclonal
Purity	Total IgG. Protein G purified in PBS pH 7.4.
MW	93 kDa
Uniprot ID	O35136
Dilution Range	WB: 1:1000, IHC-P: 1:500
Expiration Date	12 months from date of receipt.