

Product Datasheet

PAR-1 antibody (orb303270)

Description	Chicken polyclonal antibody to PAR-1
Species/Host	Gallus
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	KLH-conjugated synthetic peptide targeted against PAR1. Accession number P25116.
Form/Appearance	Liquid
Storage	Store at -20° C.
Note	For research use only
Application notes	Par 1 (Protease-activated receptor 1) stimulate phosphoinositide hydrolysis. Specific for platelets and vascular endothelial cells. Alternative names: coagulation factor II receptor, thrombin receptor.
Clonality	Polyclonal
Purity	Purified, total IgY (chicken egg yolk immunoglobulin) in 0.9% NaCl, 0.02% sodium azide.
MW	47,4 kDa
Uniprot ID	P25116
Dilution Range	ELISA: 1:500, FACS: 1:20
Expiration Date	12 months from date of receipt.