

Product Datasheet

FN antibody (orb303369)

Description

Chicken polyclonal antibody to FN

Species/Host

Gallus

Reactivity

Human

Conjugation

Biotin

Tested

WB

Applications
Immunogen

Purified, full length native protein UniProt:P02751

Form/Appearance

Liquid

Concentration

1 mg/ml

Storage

Store at 4°C; make aliquots to avoid working with a stock. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.

Note

For research use only

Application notes

Fibronectin is a glycoprotein present in a soluble dimeric form in plasma, and in a dimeric or multimeric form at the cell surface and in extracellular matrix. Fibronectin is involved in cell adhesion and migration processes including embryogenesis, wound healing, blood coagulation, host defense, and metastasis. It occurs as dimer of two 200 kDa subunits. Fibronectin has binding domains for bacterial proteins, collagens, heparin-like molecules and fibrin. Cellular fibronectin is widely distributed in the stroma of malignant tumours.

Clonality

Polyclonal

Purity

Immunogen affinity purified IgY in 0.15M sodium chloride, 0.02M sodium phosphate, 0.1% sodium azide., pH 7.2.

MW

263 kDa

Uniprot ID
P02751
Dilution Range

WB: 1: 5000

Expiration Date

12 months from date of receipt.