




Product Datasheet

TSN antibody (orb30379)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	TSN
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	26183
Uniprot ID	Q15631
NCBI	NP_001248330.1 , NP_004613.1
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescence
analysis of ...



Immunohistochemical
analysis of formalin...