


Product Datasheet

CD134 antibody (orb303967)

Description	Mouse monoclonal antibody to CD134
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-Fr, IHC-P, IP, WB
Immunogen	HTLV 1-transformed HUT-102 cells
Target	CD134
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 3-6 µg/ml. Immunohistochemistry: Recommended dilution: 5-10 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P43489
Entrez	7293
Dilution Range	Flow cytometry: Recommended dilution: 3-6 µg/ml. Immunohistochemistry: Recommended dilution: 5-10 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...