

Product Datasheet

MEPE antibody (orb30399)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target MEPE

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 58419

Uniprot ID [Q9NQ76](#)

NCBI [NP_064588.1](#), [NP_001171626.1](#), [NP_001171624.1](#), [NP_001171623.1](#), [NP_001171625.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

