



Product Datasheet

BCL2 antibody (orb304038)

Description	BCL2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human BCL2. Please contact us for the exact immunogen sequence. The peptide is available as orb374881.
Target	BCL2
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	P10415, P49950, P10417
NCBI	21, 000633, 000624
Entrez	596
Dilution Range	IHC-P: 1:50-200 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



IHC-P image of rat bone marrow tissue
US...