



## Product Datasheet

### NIS antibody (orb304039)


|                            |                                                                                                                                                     |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | NIS antibody                                                                                                                                        |
| <b>Species/Host</b>        | Rabbit                                                                                                                                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                        |
| <b>Tested Applications</b> | ELISA, ICC, IF, IHC-P                                                                                                                               |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human NIS. Please contact us for the exact immunogen sequence. The peptide is available as orb374882. |
| <b>Target</b>              | NIS                                                                                                                                                 |
| <b>Form/Appearance</b>     | 10 mM PBS, 0.02% sodium azide                                                                                                                       |
| <b>Concentration</b>       | - 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml                                                                                      |
| <b>Storage</b>             | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.                                      |
| <b>Note</b>                | For research use only                                                                                                                               |
| <b>Isotype</b>             | IgG                                                                                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                                                                          |
| <b>Purity</b>              | Polyclonal antibodies are purified by peptide affinity chromatography                                                                               |
| <b>Uniprot ID</b>          | <a href="#">Q92911</a> , <a href="#">Q63008</a> , <a href="#">Q99PN0</a>                                                                            |
| <b>NCBI</b>                | <a href="#">11</a> , <a href="#">21</a> , <a href="#">000444</a> , <a href="#">000453</a>                                                           |
| <b>Entrez</b>              | <a href="#">6528</a>                                                                                                                                |
| <b>Dilution Range</b>      | IF/ICC: 1:100-500, IHC-P: 1:100-500                                                                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                     |



IHC-P image of human thyroid cancer tiss...



Immunohistochemical staining of paraffin...



IHC-P image of mouse ovary vessel tissue...