



Product Datasheet

Collagen III antibody (orb304041)


Description	Collagen III antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Gallus, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Collagen III. Please contact us for the exact immunogen sequence. The peptide is available as orb374884.
Target	Collagen III
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	P08121, P13941, P02461
NCBI	775628, 003656
Entrez	1281
Dilution Range	IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



IHC-P image of human thyroid cancer tiss...



IHC-P staining of rat pancreas tissue us...