

Product Datasheet

SH-4-54 (orb304127)

Description	SH-4-54 is a potent STAT inhibitor with KD of 300 nM and 464 nM for STAT3 and STAT5, respectively.
Note	For research use only
Formula	$C_{29}H_{27}F_5N_2O_5S$
MW	610.59
CAS Number	[1456632-40-8]
SMILES	<chem>C(O)(=O)C1=CC=C(N(CC2=CC=C(C3CCCCC3)C=C2)C(CN(C)S(C2=C(F)C(F)=C(F)C(F)=C2F)(=O)=O)=O)C=C1</chem>
Expiration Date	12 months from date of receipt.