



Product Datasheet

VSGP/F-spondin antibody (orb30418)

Description	Rabbit polyclonal to VSGP is also known as MGC10724 antibody, SPON 1 antibody. It promotes the outgr...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	SPON1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	90973
Uniprot ID	Q9HCB6
NCBI	NP_006099.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...