


Product Datasheet

SR1001 (orb304211)

Description	SR1001 is a high-affinity synthetic ligand of ROR α and ROR γ ; suppresses TH17 cell differentiation ...
Note	For research use only
Formula	C ₁₅ H ₁₃ F ₆ N ₃ O ₄ S ₂
MW	477.4
CAS Number	[1335106-03-0]
SMILES	CC(NC1=NC(C)=C(S(=O)(NC2=CC=C(C(C(F)(F)F)(O)C(F)(F)F)C=C2)=O)S1)=O
Expiration Date	12 months from date of receipt.



Chemical
structure of
SR1001