


## Product Datasheet

### CIQ (orb304262)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | CIQ is a positive allosteric modulator selective for NR2C/D containing NMDA receptors (EC50 = 2.7 an... |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Formula</b>         | C <sub>26</sub> H <sub>26</sub> NO <sub>5</sub> Cl                                                      |
| <b>MW</b>              | 467.94134                                                                                               |
| <b>CAS Number</b>      | [486427-17-2]                                                                                           |
| <b>SMILES</b>          | <chem>O=C(C1=CC=CC(Cl)=C1)N2C(COC3=CC=C(OC)C=C3)C4=C(C=C(OC)C(OC)=C4)CC2</chem>                         |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                         |



Chemical  
structure of  
CIQ