



Product Datasheet

RNH1 antibody (orb30444)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	RNH1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49973
Uniprot ID	P13489
NCBI	NP_002930.2 , NP_976319.1 , NP_976323.1 , NP_976317.1 , NP_976321.1 , NP_976320.1 , NP_976318.1 , NP_976322.1
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescence
analysis of ...