

Product Datasheet

EYA1 antibody (orb30453)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	EYA1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	64593
Uniprot ID	Q99502
NCBI	NP_742056.1 , NP_000494.2 , NP_742055.1 , NP_742057.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

