



Product Datasheet

PTP1B (phospho-S50) antibody (orb304530)

Description	Rabbit polyclonal antibody to PTP1B (phospho-S50)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Primate, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S50 of human PTP1B protein. The exact sequence is proprietary.
Target	PTPN1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P18031, P20417, P35821
Entrez	19246, 24697, 5770
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Western blot analysis of PTP1B (Phospho-S50) antibody



Immunohistochemical analysis of PTP1B (Phospho-S50) antibody