


Product Datasheet

Ryanodine Receptor 2 (phospho-S2808) antibody (orb304574)

Description

Rabbit polyclonal antibody to Ryanodine



Western blot analysis of Ryanodine Recep...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rabbit, Rat

Conjugation

Unconjugated

Tested Applications

IF, IH, WB

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S2808 of human Ryanodine Receptor 2 protein. The exact sequence is proprietary.

Target

RYP2

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[B0LPN4](#), [E9Q401](#), [Q92736](#)

Entrez


[6262](#), [20191](#), [689560](#)

Dilution Range


WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500

Expiration Date

12 months from date of receipt.



Immunohistochemical analysis of Ryanodin...



Immunofluorescent analysis of Ryanodine

...