

Product Datasheet

alpha synuclein (phospho-S129) antibody (orb304600)

Description

Rabbit polyclonal antibody to alpha synuclein

Species/Host

Rabbit

Reactivity

Human, Mouse, Primate, Rat, Sheep

Conjugation

Unconjugated

Tested

Applications

IF, IH, WB

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S129 of human Alpha-synuclein protein. The exact sequence is proprietary.

Target

SNCA

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[P37377](#), [O55042](#), [P37840](#)

Entrez


[29219](#), [6622](#), [20617](#)

Dilution Range


WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500

Expiration Date


12 months from date of receipt.



Western blot analysis of Alpha-synuclein...



Immunohistochemical analysis of Alpha-sy...



Immunofluorescent analysis of Alpha-synu...